

Synthetic Cannabinoids (K2) - Single Drug Urine Test Panel

Catalogue No. See Box label

For in vitro diagnostic use

The Synthetic Cannabinoids Urine Test Panel is a competitive binding, lateral flow immunochromatographic assay for qualitative and simultaneous detection of Synthetic Cannabinoids in human urine at a specified cutoff level for use in employment and insurance testing.

The test provides only preliminary test results. A more specific alternative analytical method should be used in order to obtain a confirmed result. Gas Chromatography/Mass Spectrometry (GC/MS) and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS-MS) are preferred confirmatory methods.

The test is not intended to distinguish between prescription drug or illicit drug use.

Professional judgment should be exercised with any drug test result, particularly when the preliminary result is positive.

WHAT IS SYNTHETIC CANNABINOIDS URINE TEST PANEL?

The Synthetic Cannabinoids Urine Test Panel is an immunochromatographic assay for the qualitative determination of Synthetic Cannabinoids in human urine.

WHAT IS THE CUT-OFF VALUE AND APPROXIMATE DETECTION TIME?

Drug (Identifier)	Calibrator	Cut-off level	Minimum detection time	Maximum detection time
Synthetic Cannabinoids (K2)	JWH-018 Pentanoic Acid JWH-073 Butanoic Acid	50 ng/mL 50 ng/mL	8-12 hours	Up to 5+ days

WARNINGS AND PRECAUTIONS

- Not to be used for clinical diagnosis.
- This kit is for external use only. Do not swallow.
- Discard after first use. The test cannot be used more than once.
- Do not use test kit beyond expiry date.
- Do not use the kit if the pouch is punctured or not sealed.
- Keep out of the reach of children.
- Do not read after 5 minutes.

CONTENT OF THE KIT

- Test devices, one test in one pouch. One pouch contains a test and a desiccant. The desiccant is for storage purposes only, and is not used in the test procedures.
- Leaflet with instructions for use.

MATERIAL REQUIRED BUT NOT PROVIDED

- Urine collection cup
- Timer or clock

STORAGE AND STABILITY

Store at 4°C-30°C (40°F-86°F) in the sealed pouch up to the expiration date. Keep away from direct sunlight, moisture and heat. DO NOT FREEZE.

SPECIMEN COLLECTION

WHEN TO COLLECT URINE FOR THE TEST?

The minimum detection time is 8-12 hours, urine samples may be collected 8-12 hours after the suspected drug use.


HOW TO COLLECT URINE?

- Urinate directly into the urine collection cup. Urine samples may be refrigerated at 2°C-8°C (36°F-47°F) and stored up to forty-eight hours. For longer storage, freeze the samples at -20°C (-4°F) or below.
- Bring frozen or refrigerated samples to room temperature before testing. Previously frozen or refrigerated samples should be well-mixed before analysis. Cloudy specimens should be centrifuged before analysis.
- Use only clear aliquots for testing.

TEST PROCEDURE

Test should be in room temperature 18°C-30°C (65°F-86°F)

- Open the sealed pouch by tearing along the notch. Remove the test device from the pouch.
- Hold one side of the device with one hand. Use the other hand to pull out the cap and expose the absorbent end.
- Immerse the absorbent end into the urine sample for approximately 10 seconds. **Make sure that the urine level is not above the marked line printed on the front of the device.**
- Re-cap the device and lay it flat on a clean, dry, non-absorbent surface.
- Read the drug test results at 5 minutes. **Do not read results after 5 minutes.**



Note: Results after more than 5 minutes may be not accurate and should not be read.

READING THE RESULTS

Negative (-)

A colored band is visible in each control region and the appropriate test region. It indicates that the concentration of the corresponding drug of that specific test zone is zero or below the detection limit of the test.


Preliminary positive (+)

A colored band is visible in each control region. No color band appears in the appropriate test region. It indicates a preliminary positive result for the corresponding drug of that specific test zone.

Invalid

If a colored band is not visible in each of the control region or a color band is only visible in each of the test region, the test is invalid. Another test should be run to re-evaluate the specimen. If test still fails, please contact the distributor with the lot number.

Note: There is no meaning attributed to line color intensity or width.



A preliminary positive test result does not always mean a person took drugs and a negative test result does not always mean a person did not take drugs. There are a number of factors that influence the reliability of drug tests.

IMPORTANT: The result you obtained is called preliminary for a reason. The sample should be tested by a laboratory in order to determine if a drug is actually present. Send any sample which does not give a negative result to a laboratory for further testing.

What Is A False Positive Test?

The definition of a false positive test would be an instance where a substance is identified incorrectly by the Synthetic Cannabinoids Urine Test Panel. The most common causes of a false positive test are cross reactants. Certain foods and medicines, diet plan drugs and nutritional supplements may cause a false positive test result with this product.

What Is A False Negative Test?

The definition of a false negative test is that the initial drug is present but isn't detected by the Synthetic Cannabinoids Urine Test Panel. If the sample is diluted, or the sample is adulterated that may cause false negative result.

TEST LIMITATIONS

- This test has been developed for testing urine samples only. No other fluids have been evaluated. DO NOT use this device to test anything but urine.
- Adulterated urine samples may produce erroneous results. Strong oxidizing agents such as bleach (hypochlorite) can oxidize drug analytes. If a sample is suspected of being adulterated, obtain a new sample.
- This test is a qualitative screening assay. It is not designed to determine the quantitative concentration of drugs or the level of intoxication.

Note: The test provides only preliminary test results. A more specific alternative chemical method should be used in order to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS-MS) are preferred confirmatory methods. Professional judgment should be exercised with any drug test result, particularly when the preliminary result is positive.

SUMMARY

Synthetic cannabinoids are psychoactive designer drugs derived of natural herbs sprayed with synthetic chemicals that, when consumed, allegedly mimic the effects of cannabis, they are best known by the brand names K2 and Spice. Synthetic cannabinoids act on the body in a similar way to cannabinoids naturally found in cannabis, such as THC. Although synthetic cannabinoids do not produce positive results in drug tests for cannabis, it is possible to detect the metabolites in human urine.

PRINCIPLE

The Synthetic Cannabinoids Urine Test Panel is a competitive immunoassay that is used to screen for the presence of drugs in urine. It is chromatographic absorbent device in which drugs in a sample competitively combined to a limited number of drug monoclonal antibody (mouse) conjugate binding sites.

When the absorbent end is immersed into urine specimen, the urine is absorbed into the device by capillary action, mixes with the respective drug monoclonal antibody conjugate, and flows across the pre-coated membrane. When sample drug levels are zero or below the target cutoff (the detection sensitivity of the test), respective drug monoclonal antibody conjugate binds to the respective drug-protein (duck egg) conjugate immobilized in the Test Region (T) of the device. This produces a colored Test line that, regardless of its intensity, indicates a negative result.

When sample drug levels are at or above the target cutoff, the free drug in the sample binds to the respective drug monoclonal antibody conjugate preventing the respective drug monoclonal antibody conjugate from binding to the respective drug-protein conjugate immobilized in the Test Region (T) of the device. This prevents the development of a distinct colored band in the test region, indicating a potentially positive result.

To serve as a procedure control, a colored line will appear at the Control Region (C), where the Goat anti mouse IgG polyclonal antibody immobilized in, if the test has been performed properly.

QUALITY CONTROL

Users should follow the appropriate federal, state, and local guidelines concerning the frequency of assaying external quality control materials.

Even though there is an internal procedural control line in the test device in the Control Region, the use of external controls is strongly recommended as good laboratory testing practice to confirm the test procedure and to verify proper test performance. Positive and negative controls should give the expected results. When testing the positive and negative controls, the same assay procedure should be adopted. External Control (positive and negative) should be run with each new lot of test received, each new shipment, each new operator and monthly to determine that tests are working properly. This will ensure that the end user has clear understanding of when to perform quality control testing.

PERFORMANCE CHARACTERISTICS

Accuracy

Eighty clinical urine specimens were analyzed by GC-MS and by the Synthetic Cannabinoids Urine Test Panel. Each test was read by three viewers. Samples were divided by concentration into five categories: drug-free, less than half the cutoff, near cutoff negative, near cutoff positive, and high positive. Results were as follows:

Drug test	Result	Drug -free	Less than half the cutoff concentration by GC/MS analysis	Near Cutoff Negative (Between 50% below the cutoff and the cutoff concentration)	Near Cutoff Positive (Between the cutoff and 50% above the cutoff concentration)	High Positive (greater than 50% above the cutoff concentration)	%Agreement with GC/MS (95%CI)	
								Viewer
K2	A	-	10	12	17	0	0	97.5% (82% - 100%)
	Viewer	+	0	0	0	17	22	97.5% (82% - 100%)
	B	-	10	12	18	1	0	100% (84.5% - 100%)
Viewer	+	0	0	0	15	22	92.5% (77% - 100%)	
	C	-	10	12	18	3	0	100% (84.5% - 100%)

Precision and Sensitivity

To investigate the precision and sensitivity, each drug sample was analyzed at the following concentrations: cutoff - 100%, cutoff - 75%, cutoff - 50%, cutoff - 25%, cutoff, cutoff +25%, cutoff + 50%, cutoff + 75% and the cutoff + 100%. All concentrations were confirmed with GC-MS. The study was performed 2 runs/day and lasted 25 days using three different lots of the Synthetic Cannabinoids Urine Test Panel. Totally 3 operators participated in the study of the Synthetic Cannabinoids Urine Test Panel. Each of the 3 operators tests 2 aliquots at each concentration for each lot per day (2 runs/day), for a total of 50 determinations per concentration per lot of the Synthetic Cannabinoids Urine Test Panel.

Drug test	Approximate concentration of sample (ng/mL)	Number of determinations per lot	Results Negative/ Positive		
			Lot 1	Lot 2	Lot 3
K2 JWH-018 Pentanoic Acid	0	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	25.0	50	50/0	50/0	50/0
	37.5	50	50/0	50/0	50/0
	50.0	50	5/45	6/44	5/45
	62.5	50	0/50	0/50	0/50
	75.0	50	0/50	0/50	0/50
	87.5	50	0/50	0/50	0/50
K2 JWH-073 Butanoic Acid	0	50	50/0	50/0	50/0
	12.5	50	50/0	50/0	50/0
	25.0	50	50/0	50/0	50/0
	37.5	50	50/0	50/0	50/0
	50.0	50	5/45	6/44	5/45
	62.5	50	0/50	0/50	0/50
	75.0	50	0/50	0/50	0/50
	87.5	50	0/50	0/50	0/50
100.0	50	0/50	0/50	0/50	

Specificity and Cross Reactivity

To test the specificity of the test, the test device was used to test various drugs, drug metabolites and other components of the same class that are likely to be present in urine. All the components were added to drug-free normal human urine. The following structurally related compounds produced positive results with the test when tested at levels equal to or greater than the concentrations listed below.

Items	Concentration (ng/mL)
Synthetic Cannabinoids (K2)	
JWH-018 Pentanoic Acid	50
JWH-073 Butanoic Acid	50
JWH-018 N-4-hydroxypentyl	2,000
JWH-018 (Spice Cannabinoid)	1,000
JWH-018 4-Hydroxypentyl metabolite-D5 (indole-D5)	1,000
JWH-073 (Spice Cannabinoid)	2,000
JWH-073 3-Hydroxybutyl metabolite	1,000
JWH-073 3-Hydroxybutyl metabolite-D5 (indole-D5)	1,000
JWH-019 6-hydroxypentyl	1,000
JWH-122 N-4-hydroxypentyl	2,000
JWH-210 5-Hydroxypentyl metabolite	5,000
AM2201 4-Hydroxypentyl metabolite	1,000

Effect of Urinary Specific Gravity

12 urine samples with density ranges (1.005-1.025) were collected and spiked with Synthetic Cannabinoids at 25% below and 25% above cutoff levels. Each sample was tested by three batches of the Synthetic Cannabinoids Urine Test Panel. Three laboratory assistants read the result per batch of the Synthetic Cannabinoids Urine Test Panel. The results demonstrate that varying ranges of urinary specific gravity do not affect the test result

Effect of Urinary PH

The pH of an aliquot of negative urine pool was adjusted to a pH range of 4 to 9 in 1 pH unit increments and spiked with Synthetic Cannabinoids at 25% below and 25% above cutoff levels. Each sample was tested by three batches of the Synthetic Cannabinoids Urine Test Panel. Three laboratory assistants read the result per batch of the Synthetic Cannabinoids Urine Test Panel. The result demonstrates that varying range of pH do not interfere with the performance of the test.

Interfering Substances

Clinical urine samples may contain substances that could potentially interfere with the test. The following compounds were added to drug-free urine, urine with a drug concentration 25% below the cutoff, and urine with a drug concentration 25% above the cutoff for Synthetic Cannabinoids (K2). All potential interferences were added at a concentration of 100 µg/mL. None of the urine samples tested showed any deviation from the expected results.

Acetaminophen	Estroven	Nitroglycerin
Acid Reducer	Fenofibrate	Noscapine
Acyclovir	Fluoxetine HCl	Olanzapine
Advil	Fluoxamine	Omeprazole
Aleve	Fuel	Paliperidone
Amiodarone HCl	Gabapentin	Papaverine
Amlodipine Mesylate	Glibenclamide	Paroxetine
Amoxicillin	Gliclazide	Penfluridol
Ampicillin	Glipizide	Penicillin V Potassium
Aripiprazole	Glucose	Pethidine HCl
Aspirin	Haloperidol	Pioglitazone HCl
Atorvastatin Calcium	Heartburn Relief	Piracetam
Atropine	Hydrochlorothiazide	Pravastatin Sodium
Benadryl	I Caps	Prednisone Acetate
Captopril	Isosorbide Esters	Propranolol HCl
Carbamazepine	Ketoconazole	Propylthiouracil
Cefaclor	Lamotrigine	Pseudoephedrine HCl
Cephalexin	Lansoprazole	Quetiapine
Cephadrine	Levofloxacin	Risperidone
Ciprofloxacin HCl	Levonorgestrel	Ritampicin
Clarithromycin	Levothyroxine Sodium	Sertraline HCl
Clopidogrel Bisulfate	Lidocaine HCl	Sildenafil Citrate
Clozapine	Lisinopril	Simvastatin
Cortisone	Lithium Carbonate	Spiroonolactone
CVS	Loratadine	Tetracycline
Dextromethorphan HBr	Magnesium	Topiramate
Diclofenac	Mega-T Plus	Trazodone HCl
Digoxin	Metoprolol Tartrate	Triamterene
Diphenoxylate HCl	Mifepristone	Valproate
Dirithromycin	Mirtazapine	Venlafaxine HCl
Domperidone	Montelukast Sodium	Vitamin B1
Duloxetine	Mosapride	Vitamin B2
Enalapril Maleate	Nifedipine	Vitamin C
Epinephrine HCl	Nikethamide	Zencore Plus2
Esomeprazole Magnesium	Nimodipine	
Estrogen		

BIBLIOGRAPHY OF SUGGESTED READING

Baselt, R.C. Disposition of Toxic Drugs and Chemicals in Man. Biomedical Publications, Davis, CA, 1982.
Ellenhorn, M.J. and Barceloux, D. G Medical Toxicology. Elsevier Science Publishing Company, Inc., New York, 1988
Gilman, A. G., and Goodman, L. S. The Pharmacological Fluids, in Martin WR(ed): Drug Addiction I, New York, Spring – Verlag, 1977.
Harvey, R.A., Champe, P.C. Lippincotts Illustrated Reviews. Pharmacology. 91-95, 1992.
Hawwks RL, CN Chiang. Urine Testing for drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monography 73, 1986
Hofmann F.E., A Handbook on Drug and Alcohol Abuse: The Biomedical Aspects, New York, Oxford University Press, 1983.
McBay, A. J. Clin. Chem. 33,33B-40B, 1987.

INDEX OF SYMBOLS



Keep away from sunlight



Store between 4°C - 30°C (40°F - 86°F)



Keep dry



Do not re-use



In vitro diagnostic use

Distributed by:
Alere Toxicology Services
Portsmouth, VA 23704 USA