

Professional accuracy in a convenient, easy-to-use device.

# AlcoMate<sup>®</sup> Premium

Breath Alcohol Tester (Model AL7000)


- » No Calibration (PRISM<sup>™</sup> Sensor Technology)
- » FDA 510(k) Cleared
- » DOT/NHTSA Approved
- » Meets US Coast Guard Requirements
- » ISO 9002 Certified



Substance abuse testing with more **substance**.

## Features & Benefits

- » No calibration service required (PRISM™ Sensor Technology)
- » DOT approved for law enforcement (accuracy +/-0.01 at 0.10% BAC)
- » Easy one-button operation
- » Easy-replace PRISM sensor panel and sensor modules
- » Back-lit 4-digit (0.000%) LCD display
- » Active pressure sensor / deep-lung air Sampling
- » Cumulative test counter / audible alerts / low battery indicator
- » One-way-flow mouthpieces / internal mouthpiece storage slot
- » Self-diagnostic system / auto power off



## AlcoMate® Components

- |                   |                           |
|-------------------|---------------------------|
| 1 Power button    | 5 PRISM™ (Sensor Module)  |
| 2 Digital display | 6 Battery cover           |
| 3 Mouthpiece      | 7 Easy-replace panel      |
| 4 Mouthpiece slot | 8 Mouthpiece storage slot |



NOTE: The only legally valid measure of %BAC (percent Blood Alcohol Concentration) is obtained by blood testing or other evidential method of measurement. This device is for screening purposes only and should be used as a monitoring guide. This device must be used strictly in accordance with the instructions provided with the product.

Product	Description	Alere Item #
<b>AlcoMate® Premium Kit DOT Approved</b>	Kit includes: Breathalyzer, User's Manual, 2 x AA Batteries, 55 x mouthpieces, 1 pre-installed Sensor Module (PRISM™), 1 additional Sensor Module, hard plastic carrying case, car adapter, handstrap	AL 7000F-Kit
<b>AlcoMate® Sensor Module</b>	1 Pre-calibrated replacement Sensor Module for AlcoMate Premium device	PRISM
<b>AlcoMate® Mouthpieces</b>	50-Pack of standard mouthpieces for AlcoMate Premium device	ATM-Mouthpiece