




Fast, Easy 2-minute Screening Test
For the Detection of Alcohol

by



TM

ISO 13485 Registered


**FDA
CLEARED**

**CLIA
WAIVED**

The first saliva alcohol screening test on the market, Alco-Screen is a simple, cost-effective solution to detect alcohol use. Developed in 1984, you can trust Alco-Screen's quality and rely on Chematics as your partner in your alcohol testing program.



MADE IN USA



ADVANTAGES

Economical

Convenient
use anywhere, anytime

Quick
2-minute test

FDA Cleared

CLIA Waived

Over the Counter use

USES

Zero-tolerance testing

Workplace testing


**Alcohol treatment
programs**

**Criminal justice
testing**


**Emergency room
admission testing**



PROCEDURE



1. Apply saliva directly to the test pad.



2. At two minutes, compare test pad with the color chart printed on back panel of the foil package. Within 30 seconds choose the closest color match to estimate the alcohol concentration.

Contact your dealer for more information

MADE IN USA

55621 Rev 1

