

Distributed By:
Blue Grass Drug Screen, Inc.
www.uscreeningsource.com
1-866-323-7336

14 panel CLIA-waived drug test cup.

ToxCup



Urine Drug Screening Device.





For in vitro diagnostic use only.

How it works.

PRODUCT TRAINING AND CERTIFICATION AVAILABLE ONLINE.

Features & Benefits

- » CLIA-waived
- » Integrated test cup targets a wide range of illicit and prescription drugs
- » Simple procedure: collect specimen, tilt cup and read results in five minutes
- » Buprenorphine screening combined with low-level opiate cutoff makes it ideal for clinical settings
- » Self-contained cup is ideal for sending presumptive positive specimens for confirmation
- » Long shelf-life



IMPORTANT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings. For in vitro diagnostic use only. Urine drug screen assay's provide only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) or liquid chromatography/tandem mass spectrometry (LC/MS/MS) are the preferred confirmatory methods. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

Drug Test Configuration

Item #	iCup® Dx 14
DT-14	AMP500, BAR300, BUP10, BZO300, COC150, MAMP500, MDMA500, MTD300, OPI300, OXY100, PCP25, PPX300, THC50, TCA1000