

Confront workplace substance abuse with a relevant test.

UScreen® Cup

Urine Drugs of Abuse Test

Discover expanded drug detection for present-day drugs of abuse including: Fentanyl, K2/Spice, and Tramadol.

For in vitro diagnostic use. For use in employment and insurance testing.




UScreen® Cup Urine Drugs of Abuse Test

The UScreen® Cup offers fast and convenient detection for multiple drug assays in urine, and is intended for employment and insurance testing. Numerous drug configurations are available to address the challenges with workplace substance abuse. New tests help you target: misuse of prescribed Fentanyl as well as illicit use, abuse of “K2/Spice” (JHW-073 and JWH-018 synthetic cannabinoids), and emerging concerns with Tramadol abuse.

- Screens for a wide range of drugs: AMP, BAR, BUP, BZO, COC, FTY, K2, MDMA, MET, MTD, OPI, OXY, THC, TRA*
- Simple to use, no need to tilt cup to activate test—read results in five minutes
- Self-contained cup is ideal for sending preliminary positive specimens to the lab for confirmation
- Includes specimen validity test (SVT) and temperature strip—helps to monitor specimen tampering

Easy to Use



PRODUCT PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings. For in vitro diagnostic use. For use in employment and insurance testing. The assay provides a qualitative, preliminary test result. A more specific analytical method must be used in order to obtain a confirmed result. Gas chromatography-mass spectrometry (GC-MS) or liquid chromatography-tandem mass spectrometry (LC-MS/MS) are preferred confirmatory methods. Professional judgment should be applied to any drug test result, particularly when preliminary results are positive.

Ordering Information

Item #	UScreen® Cup with Specimen Validity Test	Packaging
USSCupA-10FEN	10 Panel: AMP, BUP, COC, FTY, K2, MET, OPI2000, OXY, THC, TRA (CR, PH, OX)	25 tests per box
USSCupA-10K2	10 Panel: AMP, BAR, BUP, BZO, COC, K2, MET, OPI2000, OXY, THC (CR, PH, OX)	25 tests per box
USSCupA-13FEN	13 Panel: AMP, BAR, BUP, BZO, COC, FTY, MET, MDMA, MOP300, MTD, OXY, THC, TRA (CR, PH, OX)	25 tests per box

* DRUG TESTS: Amphetamine (AMP), Barbiturates (BAR), Benzodiazepines (BZO), Buprenorphine (BUP), Cocaine (COC), Fentanyl (FTY), Marijuana (THC), Methadone (MTD), Methamphetamine (MET), Methylenedioxymethamphetamine (MDMA), Opiates (MOP/OPI), Oxycodone (OXY), Synthetic Cannabinoids (K2), Tramadol (TRA).
SPECIMEN VALIDITY TESTS: Creatinine (CR), Potential of Hydrogen (PH), Oxidants (OX).

