


# T-Cube<sup>®</sup> Oral Fluid Drug Testing



## Oral Fluid Drug Testing Benefits


Oral Fluid Drug Testing offers a simple and convenient method to perform drug screening. 100% observation during collection allows the testing to be performed anytime and anywhere, which greatly reduces the specimen-tampered concerns. Because drugs take time to metabolize and pass through body system, oral fluid is also recognized as the best way to screen the recent drug use.



## Employment & Insurance (E&I) Exemption

**NOW.**

it's the time to expand the application scenario of Oral Fluid Drug Testing to **Non-Federal Workplace and Insurance!**



The Food and Drug Administration (FDA or Agency) is announcing a list of class II devices that the Agency has determined based on established factors to no longer require premarket notification to provide reasonable assurance of safety and effectiveness, subject to certain limitations. FDA is publishing this notice of that determination in accordance with procedures established by the 21st Century Cures Act.

**E&I exempt products can be legally used in Non-Federal Workplace and Insurance!**

\* The employment & Insurance (E&I) Exemption Notice by the Food and Drug Administration on 07/11/2017.

## E&I Exempt Oral Fluid Drugs and Cutoffs

Drug	Cutoff
Barbiturates (BAR)	20 ng/mL
Methylenedioxy-methamphetamine (MDMA)	50 ng/mL
Marijuana (THC)	40 ng/mL
Amphetamine (AMP)	50 ng/mL
Cocaine (COC)	20 ng/mL

Drug	Cutoff
Methamphetamine (MET)	50 ng/mL
Opiate (OPI)	40 ng/mL
Oxycodone (OXY)	20 ng/mL
Methadone (MTD)	30 ng/mL
Phencyclidine (PCP)	10 ng/mL

# T-Cube® Oral Fluid Drug Testing



- Observed and Tamper-resistant Collection
- E&I Exempt Testing Panels
- Saliva Saturation Indicator
- Fast Results

We are offering E&I exempt Oral Fluid Drug Testing products as below:

Item Code	Description	Complexity
ODOA-856EUO	5-Panel T-Cube: AMP, COC, MET, OPI, PCP	E&I Exempt
ODOA-966EUO	6-Panel T-Cube: AMP, BAR, COC, OPI, PCP, THC	E&I Exempt
ODOA-376EUO	7-Panel T-Cube: AMP, COC, MET, OPI, OXY, PCP, THC	E&I Exempt
ODOA-10106EUO	10-Panel T-Cube: AMP, BAR, COC, MDMA, MET, MTD, OPI, OXY, PCP, THC	E&I Exempt

Distributed by: Blue Grass Drug Screen, Inc.  
[www.ussscreeningsource.com](http://www.ussscreeningsource.com)  
 1-866-323-7336